

“The Spooner” Patella Stabilizing Strapping System Directions for Use



1. Lay brace on top of “Spooner”, felt side up and short arms on the lateral side of the brace.



2. Wrap and attach short arms above and below the lateral hinge plate and attach to the “Spooner” with the velcro hook squares. Cut excess as needed.



3. Wrap and attach long arms above and below the medial hinge plate and attach to the “Spooner” with the alligator clips, overlapping the previously attached short arms. Cut excess as needed.



1850 DIPLOMAT DRIVE #100, FARMERS BRANCH, TX 75234
OFFICE (800) 872-5488 FAX (800) 455-5488



“The Spooner” Patella Stabilizing Strapping System Directions for Use



1. Lay brace on top of “Spooner”, felt side up and short arms on the lateral side of the brace.



2. Wrap and attach short arms above and below the lateral hinge plate and attach to the “Spooner” with the velcro hook squares. Cut excess as needed.



3. Wrap and attach long arms above and below the medial hinge plate and attach to the “Spooner” with the alligator clips, overlapping the previously attached short arms. Cut excess as needed.



1850 DIPLOMAT DRIVE #100, FARMERS BRANCH, TX 75234
OFFICE (800) 872-5488 FAX (800) 455-5488

